



Information for Customers
About Latex in Zimmer Products
 (updated 7/5/2011)

Most products manufactured by Zimmer are free from natural rubber latex. If the item number is not listed below, it does not, nor has it ever, contained natural rubber latex.

Only products that either contain, or have contained in the past, natural rubber latex are listed below. If applicable, the date or lot number that the product changed to latex-free is listed as a comment. If no comment exists then the product currently contains natural rubber latex.

There are no products from the Zimmer Austin, Dental, GmbH, Spine or Trabecular Metal Technology divisions that contain natural rubber latex.

If more information is needed, please contact Zimmer, Inc., P.O. Box 708, Warsaw, IN 46580 (Attn: Regulatory Affairs), Phone: 574-267-6131.

Division	Item Number	Description	Comments
Zimmer Warsaw	00-5049-091-00	<i>Osteobond</i> [®] Vacuum Mixing System	Lot number 31255700 and higher are latex free. Earlier systems had a small rubber band attached to the tubing set.
	00-5049-061-00	Vacuum Mixing Bowl	Lot number 531799 and higher are latex free. Earlier systems had a small rubber band attached to the tubing set.
Zimmer Orthopaedic Surgical Products	00-0766-001-00	Ankle Brace with Elastic Strap, X-Small, Left	Devices manufactured later than 10/16/2000 are latex free.
	00-0766-002-00	Ankle Brace with Elastic Strap, X-Small, Right	Devices manufactured later than 10/16/2000 are latex free.
	00-0766-003-00	Ankle Brace with Elastic Strap, Small, Left	Devices manufactured later than 10/16/2000 are latex free.
	00-0766-004-00	Ankle Brace with Elastic Strap, Small, Right	Devices manufactured later than 10/16/2000 are latex free.
	00-0766-005-00	Ankle Brace with Elastic Strap, Medium, Left	Devices manufactured later than 10/16/2000 are latex free.
	00-0766-006-00	Ankle Brace with Elastic Strap, Medium, Right	Devices manufactured later than 10/16/2000 are latex free.
	00-0766-007-00	Ankle Brace with Elastic Strap, Large, Left	Devices manufactured later than 10/16/2000 are latex free.
	00-0766-008-00	Ankle Brace with Elastic Strap, Large, Right	Devices manufactured later than 10/16/2000 are latex free.

Division	Item Number	Description	Comments
Zimmer Orthopaedic Surgical Products (continued)	00-0766-009-00	Ankle Brace with Elastic Strap, X-Large, Left	Devices manufactured later than 10/16/2000 are latex free.
	00-0766-010-00	Ankle Brace with Elastic Strap, X-Large, Right	Devices manufactured later than 10/16/2000 are latex free.
	00-1799-001-00	Adjustable Wrist Brace, X-Small, Left	Devices manufactured later than 2/25/2002 are latex free.
	00-1799-002-00	Adjustable Wrist Brace, X-Small, Right	Devices manufactured later than 2/25/2002 are latex free.
	00-1799-003-00	Adjustable Wrist Brace, Small, Left	Devices manufactured later than 2/25/2002 are latex free.
	00-1799-004-00	Adjustable Wrist Brace, Small, Right	Devices manufactured later than 2/25/2002 are latex free.
	00-1799-005-00	Adjustable Wrist Brace, Medium, Left	Devices manufactured later than 2/25/2002 are latex free.
	00-1799-006-00	Adjustable Wrist Brace, Medium, Right	Devices manufactured later than 2/25/2002 are latex free.
	00-1799-007-00	Adjustable Wrist Brace, Large, Left	Devices manufactured later than 2/25/2002 are latex free.
	00-1799-008-00	Adjustable Wrist Brace, Large, Right	Devices manufactured later than 2/25/2002 are latex free.
	00-1799-009-00	Adjustable Wrist Brace, X-Large, Left	Devices manufactured later than 2/25/2002 are latex free.
	00-1799-010-00	Adjustable Wrist Brace, X-Large, Right	Devices manufactured later than 2/25/2002 are latex free.
	00-2500-000-10	400ml Compact Evacuator Kit, 1/8" PVC with Trocar	Lot number 32285900 and higher are latex free.
	00-2500-004-10	400ml Compact Evacuator Component Kit	Lot number 32286200 and higher are latex free.
	00-2500-700-10	400ml Compact Evacuator, 1/8" Drain with Bonded Suction Port	Lot number 32546300 and higher are latex free.
	00-2500-900-10	400ml Compact Evacuator, 1/8" PVC Hydrocoat with Trocar	Lot number 32360300 and higher are latex free.
	00-2505-000-10	400ml Compact Evacuator, 3/16" PVC with Trocar	Lot number 32286300 and higher are latex free.
	00-2505-700-10	400ml Compact Evacuator, 3/16" Drain with Bonded Suction Port	Lot number 32302400 and higher are latex free.
	00-2540-000-10	400ml Compact Evacuator, 1/8" Middle Perf Silicone with Trocar	Lot number 32342600 and higher are latex free.

Division	Item Number	Description	Comments
Zimmer Orthopaedic Surgical Products (continued)	00-2541-000-10	400ml Compact Evacuator, 3/16" Middle Perf Silicone with Trocar	Lot number 32308800 and higher are latex free.
	00-2545-000-10	400ml Compact Evacuator, 1/4" Middle Perf Silicone with Trocar	Lot number 32308900 and higher are latex free.
	00-2550-001-10	Infection Control Kit with 3/32" Drain	Lot number 32328200 and higher are latex free.
	00-2550-002-10	Infection Control Kit with 1/8" Drain	Lot number 32328300 and higher are latex free.
	00-2550-003-10	Infection Control Kit with 3/16" Drain	Lot number 32332000 and higher are latex free.
	00-2550-004-10	Infection Control Kit with 1/4" Drain	Lot number 32401100 and higher are latex free.
	00-2550-050-10	Infection Control Compact Kit	Lot number 32355700 and higher are latex free.
	00-2550-703-10	400ml Infection Control Evacuator, 3/16" Drain with Bonded Suction Port	Lot number 32486600 and higher are latex free.
	00-2550-704-10	400ml Infection Control Evacuator, 1/4" Drain with Bonded Suction Port	Lot number 32486700 and higher are latex free.
	00-2555-020-05	Hemovac Autotransfusion System, 1/8" Drain	Lot number 60051067 and higher are latex free.
	00-2555-030-05	Hemovac Autotransfusion System, 3/16" Soft Drain	Lot number 60057360 and higher are latex free.
	00-2555-040-05	Hemovac Autotransfusion System, 1/4" Drain	Lot number 60056540 and higher are latex free.
	00-2555-060-05	Hemovac Autotransfusion System, Without Drain	Lot number 60093196 and higher are latex free.
	00-2587-001-10	400ml Compact Urevac Kit with 1/16" Trocar	
	00-2590-000-10	400ml Compact Evacuator, 1/4" PVC with Trocar (4")	Lot number 32287500 and higher are latex free.
	00-2590-040-10	400ml Compact Evacuator, 1/4" PVC with Trocar	Lot number 32337700 and higher are latex free.
	00-2590-700-10	400ml Compact Evacuator, 1/4" Drain with Bonded Suction Port	Lot number 32887100 and higher are latex free.
	00-2590-900-10	Hydrocoat PVC Kit, 1/4" with Trocar	Lot number 32393300 and higher are latex free.
	00-2595-000-10	400ml Compact Evacuator, 3/32" PVC with Trocar	Lot number 32303700 and higher are latex free.
	00-3860-001-00	Flex-Foam Traction Strips, 2" x 42"	

Division	Item Number	Description	Comments
Zimmer Orthopaedic Surgical Products (continued)	00-3860-002-00	Flex-Foam Traction Strips, 3" x 42"	
	00-3860-003-00	Flex-Foam Traction Strips, 4" x 42"	
	00-3861-001-00	Flex-Foam Traction Strips, 2" x 25 yards	
	00-3861-002-00	Flex-Foam Traction Strips, 3" x 25 yards	
	00-3861-003-00	Flex-Foam Traction Strips, 4" x 25 yards	
	00-4444-002-00	Standard Orthopaedic Wraps, 2 x 5 yards	
	00-4444-003-00	Standard Orthopaedic Wraps, 3 x 5 yards	
	00-4444-004-00	Standard Orthopaedic Wraps, 4 x 5 yards	
	00-4444-006-00	Standard Orthopaedic Wraps, 6 x 5 yards	
	00-4444-012-00	Standard Orthopaedic Wraps, 2 x 5 yards	
	00-4444-013-00	Standard Orthopaedic Wraps, 3 x 5 yards	
	00-4444-014-00	Standard Orthopaedic Wraps, 4 x 5 yards	
	00-4444-016-00	Standard Orthopaedic Wraps, 6 x 5 yards	
	00-4444-024-00	Standard Orthopaedic Wraps, 4" x 10 yards	
	00-4444-026-00	Standard Orthopaedic Wraps, 6" x 10 yards	
	00-5150-035-00	Dual Spike Fluid Set	
	00-5160-030-00	General Surgery Tubing Set	Lot number 32947800 and higher are latex free.
	00-5190-020-00	CDIS Irrigation Fluid Set	
	00-5190-021-00	CDIS High Flow Tubing Set	
	00-5555-002-00	Premium Orthopaedic Wraps, 2" x 5-1/2 yards	
00-5555-003-00	Premium Orthopaedic Wraps, 3" x 5-1/2 yards		

Division	Item Number	Description	Comments
Zimmer Orthopaedic Surgical Products (continued)	00-5555-004-00	Premium Orthopaedic Wraps, 4" x 5-1/2 yards	
	00-5555-006-00	Premium Orthopaedic Wraps, 6" x 5-1/2 yards	
	00-5555-024-00	Premium Orthopaedic Wraps, 4" x 11 yards	
	00-5555-026-00	Premium Orthopaedic Wraps, 6" x 11 yards	
	00-5555-032-00	Sterile Premium Orthopaedic Wraps, 2" x 5-1/2 yards	
	00-5555-033-00	Sterile Premium Orthopaedic Wraps, 3" x 5-1/2 yards	
	00-5555-034-00	Sterile Premium Orthopaedic Wraps, 4" x 5-1/2 yards	
	00-5555-036-00	Sterile Premium Orthopaedic Wraps, 6" x 5-1/2 yards	
	00-5555-044-00	Sterile Premium Orthopaedic Wraps, 4" x 11 yards	
	00-5555-046-00	Sterile Premium Orthopaedic Wraps, 6" x 11 yards	
	60-2694-001-00	Esmark Bandages, Disposable 3" x 9'	
	60-2694-002-00	Esmark Bandages, Disposable 4" x 9'	
	60-2694-003-00	Esmark Bandages, Disposable 6" x 9'	
	60-2694-004-00	Esmark Bandages, Disposable 3" x 12'	
	60-2694-005-00	Esmark Bandages, Disposable 4" x 12'	
	60-2694-006-00	Esmark Bandages, Disposable 6" x 12'	
Hall/Linvatec	60-5175-036-00	Irrigation Fluid Pump Set	Latex is contained in the pressure sensing bladder scope/cannula connectors
	60-5175-038-00	Quick Connect Tubing Set	Latex is contained in the pressure sensing bladder scope/cannula connectors
	C7039	Aqua-pro Tubing Set	

Division	Item Number	Description	Comments
Hall/Linvatec (continued)	9321	Puddle-Vac	Latex is contained in the rubber band that is placed around the device during packaging.
	9350	Small Joint Tubing and Cannula Set	Lot number 58014 and higher are latex free.
	83100	Tube Set Arthro, Use with 83000 Pump	
	83100N	Tube Set Non-Sterile	
	87100	Tube Set Arthro Cassette, Use with 87000 Pump	
	83091	Tube Set Arthro, Use with Small Joint Cannula System	