



## Zimmer® Periarticular Proximal Lateral Humeral Locking Plate

### Plate Features:

- Anatomically contoured plates create a fit that requires little or no additional bending and are helpful for metaphyseal/diaphyseal reduction
- The low profile plate facilitates fixation without impinging on soft tissues
- Plates are available in a variety of sizes and lengths, from 4 to 16 holes, left and right
- Dual-compression slots accommodate periarticular screws or conventional stainless steel screws and allow bi-directional compression
- The last diaphyseal plate hole is designed to accommodate the tension device (00-4817-000-05)
- The tapered plate shaft design allows for a minimally invasive technique with submuscular passage of the plate
- Multiple locking holes in the plate head allow placement of the screws to capture fragments
- Thick-to-thin profiles make the plates autocontourable
- The anatomical shape of the head of the plate matches the shape of the proximal humerus



| Prod. No. (Left)                                    | Prod. No. (Right) | Description                            | Holes | Length |
|---|-------------------|--|-------|--------|
| <b>3.5mm Proximal Lateral Humeral Locking Plate</b> |                   |  |       |        |
| 00-2358-002-04                                      | 00-2358-001-04    | Proximal Lateral Humeral Locking Plate | 4     | 90mm   |
| 00-2358-002-06                                      | 00-2358-001-06    | Proximal Lateral Humeral Locking Plate | 6     | 114mm  |
| 00-2358-002-08                                      | 00-2358-001-08    | Proximal Lateral Humeral Locking Plate | 8     | 138mm  |
| 00-2358-002-12                                      | 00-2358-001-12    | Proximal Lateral Humeral Locking Plate | 12    | 186mm  |
| 00-2358-002-16                                      | 00-2358-001-16    | Proximal Lateral Humeral Locking Plate | 16    | 234mm  |

| Prod. No.  | Description  | Set Qty. |
|--|--|----------|
| <b>00-2358-000-05 Proximal Lateral Humeral Locking Plate Set</b> |  |          |
| 00-2358-002-04   | Proximal Lateral Humeral Locking Plate, 4 Hole 90mm Lng, Left    | 1        |
| 00-2358-002-06   | Proximal Lateral Humeral Locking Plate, 6 Hole 114mm Lng, Left   | 1        |
| 00-2358-002-08   | Proximal Lateral Humeral Locking Plate, 8 Hole 138mm Lng, Left   | 1        |
| 00-2358-002-12   | Proximal Lateral Humeral Locking Plate, 12 Hole 186mm Lng, Left  | 1        |
| 00-2358-002-16   | Proximal Lateral Humeral Locking Plate, 16 Hole 234mm Lng, Left  | 1        |
| 00-2358-001-04   | Proximal Lateral Humeral Locking Plate, 4 Hole 90mm Lng, Right   | 1        |
| 00-2358-001-06   | Proximal Lateral Humeral Locking Plate, 6 Hole 114mm Lng, Right  | 1        |
| 00-2358-001-08   | Proximal Lateral Humeral Locking Plate, 8 Hole 138mm Lng, Right  | 1        |
| 00-2358-001-12   | Proximal Lateral Humeral Locking Plate, 12 Hole 186mm Lng, Right | 1        |
| 00-2358-001-16   | Proximal Lateral Humeral Locking Plate, 16 Hole 234mm Lng, Right | 1        |

| Prod. No.   | Description                               | Set Qty. |
|---|---|----------|
| <b>00-2360-000-14 Proximal Lateral Humeral Plate Standard Jig Set</b> |   |          |
| 00-2360-092-01  | Proximal Lateral Humeral Plate Jig, Right | 1        |
| 00-2360-092-02  | Proximal Lateral Humeral Plate Jig, Left  | 1        |
| 00-2358-030-00  | Case                                      | 1        |

| Prod. No.  | Description                                      | Set Qty. |
|--|--|----------|
| <b>00-2360-000-02 3.5mm/2.7mm Periarticular Locking Instrument Set</b> |  |          |
| 00-1147-071-00   | Cleaning Stylet                                  | 1        |
| 00-1147-076-00   | Cleaning Brush                                   | 1        |
| 00-2358-040-00   | 3.5/2.7mm Locking Screw and Instrument Case      | 1        |
| 00-2360-011-01   | 3.5mm/2.7mm Locking Plate Reduction Instrument   | 2        |
| 00-2360-011-02   | 3.5mm/2.7mm Locking Plate Reduction Sleeve       | 2        |
| 00-2360-012-03   | Plate Reduction Spin Knob                        | 2        |
| 00-2360-020-20   | 2.0mm Standard Cannula                           | 4        |
| 00-2360-020-27   | 2.7mm Standard Cannula                           | 4        |
| 00-2360-021-16   | 1.6mm Standard Cannula                           | 4        |
| 00-2360-033-16   | 1.6mm Standard Drill Tip Guide Wire (5 per box)  | 2 bx     |
| 00-2360-040-35   | 3.5mm/2.7mm Locking Screw Standard Depth Gauge   | 1        |
| 00-2360-041-35   | 3.5mm Cannulated Locking Screw Depth Gauge       | 1        |
| 00-2360-053-27   | 2.7mm Locking Screw Tap                          | 2        |
| 00-2360-053-35   | 3.5mm Locking Screw Tap                          | 2        |
| 00-2360-054-35   | 3.5mm Cannulated Locking Screw Tap               | 2        |
| 00-2360-065-00   | 2.5mm Screwdriver Stop Ring                      | 1        |
| 00-2360-065-25   | 2.5mm Hex Standard Screwdriver                   | 2        |
| 00-2360-066-25   | 2.5mm Hex Standard Cannulated Screwdriver        | 2        |
| 00-2360-071-27   | 2.7mm Standard Cannulated Drill                  | 2        |
| 00-2360-080-00   | Torque Limiting Attachment                       | 1        |
| 00-2360-086-00   | Modular Handle                                   | 1        |
| 00-2360-088-00   | 3.5mm/2.7mm Cannula Insert                       | 1        |
| 00-2360-175-20   | 2.0mm Standard Drill                             | 2        |
| 00-2360-205-27   | 2.7mm Standard Drill                             | 2        |
| 00-4812-045-00   | Large Hex Screwdriver                            | 1        |
| <b>Also Available:</b>   |  |          |
| 00-2360-050-25   | 2.5mm Hex Forward Captive Screwdriver            | 0        |
| 00-2360-051-25   | 2.5mm Hex Cannulated Forward Captive Screwdriver | 0        |

**Designed for use with:**

The 3.5mm screws in the *Zimmer* 3.5mm/2.7mm Periarticular Locking Screw Set (00-2359-000-02)

The 3.5mm screws in the *Zimmer* 3.5mm/2.7mm Periarticular Locking Instrument Set (00-2360-000-02)

*Zimmer* Periarticular Proximal Lateral Humeral Locking Plate Standard Jig Set (00-2360-000-14)

**Compatible with:**

*Zimmer* Periarticular 3.5mm Self-tapping Cortical Screws

*Zimmer* Periarticular 4.0mm Cancellous Screws

*Zimmer* Small Fragment Instrument Set (00-4800-035-20)

*Zimmer* Small Fragment Screw Set (00-4800-035-01)





| Prod. No.             | Description                             | Set Qty. |                 |                              |   |
|-----------------------|---|----------|-----------------|------------------------------|---|
| <b>00-2359-000-02</b> | <b>3.5mm/2.7mm Locking Screw Set</b>    |          |                 |                              |   |
| 00-2359-030-36        | 3.5mm Cannulated Locking Screw 30mm Lng | 3        | 00-2359-040-35  | 3.5mm Locking Screw 40mm Lng | 3 |
| 00-2359-035-36        | 3.5mm Cannulated Locking Screw 35mm Lng | 3        | 00-2359-042-35  | 3.5mm Locking Screw 42mm Lng | 3 |
| 00-2359-040-36        | 3.5mm Cannulated Locking Screw 40mm Lng | 3        | 00-2359-044-35  | 3.5mm Locking Screw 44mm Lng | 3 |
| 00-2359-045-36        | 3.5mm Cannulated Locking Screw 45mm Lng | 3        | 00-2359-046-35  | 3.5mm Locking Screw 46mm Lng | 3 |
| 00-2359-050-36        | 3.5mm Cannulated Locking Screw 50mm Lng | 4        | 00-2359-048-35  | 3.5mm Locking Screw 48mm Lng | 3 |
| 00-2359-055-36        | 3.5mm Cannulated Locking Screw 55mm Lng | 4        | 00-2359-050-35  | 3.5mm Locking Screw 50mm Lng | 3 |
| 00-2359-060-36        | 3.5mm Cannulated Locking Screw 60mm Lng | 4        | 00-2359-052-35  | 3.5mm Locking Screw 52mm Lng | 2 |
| 00-2359-065-36        | 3.5mm Cannulated Locking Screw 65mm Lng | 4        | 00-2359-054-35  | 3.5mm Locking Screw 54mm Lng | 2 |
| 00-2359-070-36        | 3.5mm Cannulated Locking Screw 70mm Lng | 4        | 00-2359-056-35  | 3.5mm Locking Screw 56mm Lng | 2 |
| 00-2359-075-36        | 3.5mm Cannulated Locking Screw 75mm Lng | 3        | 00-2359-058-35  | 3.5mm Locking Screw 58mm Lng | 2 |
| 00-2359-080-36        | 3.5mm Cannulated Locking Screw 80mm Lng | 3        | 00-2359-060-35  | 3.5mm Locking Screw 60mm Lng | 2 |
| 00-2359-085-36        | 3.5mm Cannulated Locking Screw 85mm Lng | 3        | 00-2359-065-35  | 3.5mm Locking Screw 65mm Lng | 2 |
| 00-2359-090-36        | 3.5mm Cannulated Locking Screw 90mm Lng | 3        | 00-2359-070-35  | 3.5mm Locking Screw 70mm Lng | 2 |
| 00-2359-030-37        | 3.5mm Cannulated Conical Screw 30mm Lng | 2        | 00-2359-010-27* | 2.7mm Locking Screw 10mm Lng | 3 |
| 00-2359-035-37        | 3.5mm Cannulated Conical Screw 35mm Lng | 2        | 00-2359-012-27* | 2.7mm Locking Screw 12mm Lng | 3 |
| 00-2359-040-37        | 3.5mm Cannulated Conical Screw 40mm Lng | 2        | 00-2359-014-27* | 2.7mm Locking Screw 14mm Lng | 3 |
| 00-2359-045-37        | 3.5mm Cannulated Conical Screw 45mm Lng | 2        | 00-2359-016-27* | 2.7mm Locking Screw 16mm Lng | 3 |
| 00-2359-050-37        | 3.5mm Cannulated Conical Screw 50mm Lng | 2        | 00-2359-018-27* | 2.7mm Locking Screw 18mm Lng | 3 |
| 00-2359-055-37        | 3.5mm Cannulated Conical Screw 55mm Lng | 2        | 00-2359-020-27* | 2.7mm Locking Screw 20mm Lng | 3 |
| 00-2359-060-37        | 3.5mm Cannulated Conical Screw 60mm Lng | 2        | 00-2359-022-27* | 2.7mm Locking Screw 22mm Lng | 4 |
| 00-2359-065-37        | 3.5mm Cannulated Conical Screw 65mm Lng | 2        | 00-2359-024-27* | 2.7mm Locking Screw 24mm Lng | 4 |
| 00-2359-070-37        | 3.5mm Cannulated Conical Screw 70mm Lng | 2        | 00-2359-026-27* | 2.7mm Locking Screw 26mm Lng | 4 |
| 00-2359-012-35        | 3.5mm Locking Screw 12mm Lng            | 2        | 00-2359-028-27* | 2.7mm Locking Screw 28mm Lng | 4 |
| 00-2359-014-35        | 3.5mm Locking Screw 14mm Lng            | 2        | 00-2359-030-27* | 2.7mm Locking Screw 30mm Lng | 4 |
| 00-2359-016-35        | 3.5mm Locking Screw 16mm Lng            | 2        | 00-2359-032-27* | 2.7mm Locking Screw 32mm Lng | 4 |
| 00-2359-018-35        | 3.5mm Locking Screw 18mm Lng            | 2        | 00-2359-034-27* | 2.7mm Locking Screw 34mm Lng | 4 |
| 00-2359-020-35        | 3.5mm Locking Screw 20mm Lng            | 2        | 00-2359-036-27* | 2.7mm Locking Screw 36mm Lng | 4 |
| 00-2359-022-35        | 3.5mm Locking Screw 22mm Lng            | 2        | 00-2359-040-27* | 2.7mm Locking Screw 40mm Lng | 4 |
| 00-2359-024-35        | 3.5mm Locking Screw 24mm Lng            | 2        | 00-2359-042-27* | 2.7mm Locking Screw 42mm Lng | 4 |
| 00-2359-026-35        | 3.5mm Locking Screw 26mm Lng            | 2        | 00-2359-044-27* | 2.7mm Locking Screw 44mm Lng | 4 |
| 00-2359-028-35        | 3.5mm Locking Screw 28mm Lng            | 2        | 00-2359-046-27* | 2.7mm Locking Screw 46mm Lng | 4 |
| 00-2359-030-35        | 3.5mm Locking Screw 30mm Lng            | 2        | 00-2359-048-27* | 2.7mm Locking Screw 48mm Lng | 4 |
| 00-2359-032-35        | 3.5mm Locking Screw 32mm Lng            | 3        | 00-2359-050-27* | 2.7mm Locking Screw 50mm Lng | 4 |
| 00-2359-034-35        | 3.5mm Locking Screw 34mm Lng            | 3        | 00-2359-052-27* | 2.7mm Locking Screw 52mm Lng | 3 |
| 00-2359-036-35        | 3.5mm Locking Screw 36mm Lng            | 3        | 00-2359-054-27* | 2.7mm Locking Screw 54mm Lng | 3 |
| 00-2359-038-35        | 3.5mm Locking Screw 38mm Lng            | 3        | 00-2359-056-27* | 2.7mm Locking Screw 56mm Lng | 3 |
|                       |   |          | 00-2359-058-27* | 2.7mm Locking Screw 58mm Lng | 3 |
|                       |   |          | 00-2359-060-27* | 2.7mm Locking Screw 60mm Lng | 3 |

\* The proximal lateral humeral locking plate is designed for use with 3.5mm screws, and is not intended for use with the 2.7mm screw.



+H124972347108001/\$050101H072

97-2347-108-00 2.5MM Printed in USA ©2005, 2007 Zimmer, Inc.



Contact your Zimmer representative or  
visit us at [www.zimmer.com](http://www.zimmer.com)



**zimmer** | trauma  
*Confidence in your hands®*